

Table S4 Sequences for universal primers

Common primers used for attaching barcodes and sequencing adaptors. The first two primers are used in the first round of pre-sequencing PCR. The last two primers are used in the second round of pre-sequencing PCR.

Primer/purpose	Sequence
Forward common (followed by barcode and site-specific forward sequence)	CGACGCTCTCCGATCT
Reverse common (follows site-specific reverse sequence)	CATACGAGCTCTCCGATCT
Forward primer to attach sequencing adaptor	AATGATACGGCGACCACCGAGATCTACACTCTTCCCTACAG ACGCTCTCCGATCT
Reverse primer to attach sequencing adaptor	CAAGCAGAAGACGGCATAACGAGCTCTCCGATCT